

HUMAN HEALTH

ENVIRONMENTAL HEALTH

ChemBioOffice Version 14.0

Hardware & Software Guide

Last Updated: April 18, 2014



Table of Contents

1	Introduction.....	3
2	Hardware Requirements	3
3	Software Requirements.....	3
4	Technical Support	4

1 Introduction

This document is intended for general use with PerkinElmer Desktop 14.0 customers. This document does not cover details about the remote database sharing and three tier version as they are not supported for this release.

2 Hardware Requirements

Hardware Component	Minimum Specification	Recommended Specification
Memory (RAM)	For 32-bit OS 1 GB For 64-bit OS 2 GB	For 32-bit OS 1 GB For 64-bit OS 2 GB
Processor	1.6 GHz or higher	3.0 GHz or higher dual core processor
Free Space Required on Hard Drive	300 MB - 1 GB (1 GB is for installing the full ChemBioOffice Ultra). This is for installation only and does not include any creation/expansion of any databases.	10 GB
Screen resolution	1024 x 768 or higher	1024 x 768 or higher

3 Software Requirements

Software	Windows	Macintosh
Operating System	<ul style="list-style-type: none"> Windows 7 Professional and Ultimate (32 bit) Windows 7 Professional and Ultimate (64 bit) Windows 8.1 (32 bit) Windows 8.1 (64 bit) 	<ul style="list-style-type: none"> Mac OS X 10.9.x
Microsoft Office (Standard, Professional, or Enterprise)	<ul style="list-style-type: none"> Microsoft Office 2010 Microsoft Office 2013 	<ul style="list-style-type: none"> Mac OS supports only Microsoft Office 2011
.Net Framework	<ul style="list-style-type: none"> .Net Framework 4.0 <p>Important Notes:</p> <ul style="list-style-type: none"> .Net Framework 3.5 is part of the installer and should not be uninstalled. .Net Framework 3.5 and 4.0 can operate successfully side-by-side. 	N/A
Browser plug-in/ ActiveX controls:	<ul style="list-style-type: none"> Microsoft Internet Explorer 9.0 (32 and 64 bit) Microsoft Internet Explorer 10.x (32 and 64 bit) Microsoft Internet Explorer 11.x (32 and 64 bit) Mozilla Firefox 27 or higher Google Chrome 32 or higher 	<ul style="list-style-type: none"> Mac Safari 7.0 (64 bit) Mozilla Firefox 27 or higher (64 bit) Google Chrome 32 or higher

4 Technical Support

To access FAQs, KB articles, support forums and more, please visit:

<http://www.perkinelmer.com/informatics/support/>

For support contact options, please visit:

<http://www.perkinelmer.com/informatics/support/contact>

PerkinElmer Inc. | For the Better
HUMAN HEALTH | ENVIRONMENTAL HEALTH
940 Winter Street | Waltham, MA 02451