

HUMAN HEALTH

ENVIRONMENTAL HEALTH

ChemOffice Version 15.1

Hardware & Software Guide

Last Updated: March 2, 2016



Table of Contents

1	Introduction.....	3
2	Hardware Requirements	3
3	Software Requirements.....	3
4	Technical Support	4

1 Introduction

This document is intended for general use with PerkinElmer Desktop 15.1 customers.

2 Hardware Requirements

Hardware Component	Minimum Specification	Recommended Specification
Memory (RAM)	For 32-bit OS 1 GB For 64-bit OS 2 GB	For 32-bit OS 1 GB For 64-bit OS 2 GB
Processor	1.6 GHz or higher	3.0 GHz or higher dual core processor
Free Space Required on Hard Drive	300 MB - 1 GB (1 GB is for installing the full ChemOffice Professional). This is for installation only and does not include any creation/expansion of any databases.	10 GB
Screen resolution	1024 x 768 or higher	1024 x 768 or higher

3 Software Requirements

Software	Windows	Macintosh (ChemDraw Only)
Operating System	<ul style="list-style-type: none"> Windows 7 Professional and Ultimate (32 bit) Windows 7 Professional and Ultimate (64 bit) Windows 8.1 (32 bit) Windows 8.1 (64 bit) Windows 10 (32 bit) Windows 10 (64 bit) 	<ul style="list-style-type: none"> Mac OS X 10.10 and 10.11
Microsoft Office (Standard, Professional, or Enterprise) NOTE: For 32 bit only.	<ul style="list-style-type: none"> Microsoft Office 2010 Microsoft Office 2013 Microsoft Office 365 	<ul style="list-style-type: none"> Microsoft Office 2011 Microsoft Office 2016 Microsoft Office 365
.Net Framework	<ul style="list-style-type: none"> .Net Framework 4.0 and 4.5 <p>Important Notes:</p> <ul style="list-style-type: none"> .Net Framework 3.5 is part of the installer and should not be uninstalled. .Net Framework 3.5 and 4.0 can operate successfully side-by-side. 	N/A
Browser plug-in/ActiveX controls:	<ul style="list-style-type: none"> Microsoft Internet Explorer 9.x (32 and 64 bit) Microsoft Internet Explorer 10.x (32 and 64 bit) Microsoft Internet Explorer 11.x (32 and 64 bit) Mozilla Firefox 27 or higher <p>NOTE: Google Chrome 42.x or later does not support plug-ins on Windows.</p>	<ul style="list-style-type: none"> Mac Safari 9 Mozilla Firefox 28 or higher (64 bit) Google Chrome 32 through 38 <p>NOTE: Google Chrome 39.x does not support plug-ins on Mac OS.</p>

4 Technical Support

To access FAQs, KB articles, support forums and more, please visit:

<http://www.perkinelmer.com/informatics/support/>

For support contact options, please visit:

<http://www.perkinelmer.com/informatics/support/contact>

PerkinElmer Inc. | For the Better
HUMAN HEALTH | ENVIRONMENTAL HEALTH
940 Winter Street | Waltham, MA 02451